# Zirlux<sup>®</sup> Wax

Wax Discs

Discos de cera

Cire disque

Wachs-Scheiben

Cera disco

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Ø 98,5 mm

H: 20 mm

DE Healthcare Ltd  
Unit 9, Kingsthorpe Business Centre,  
Studland Road, Kingsthorpe,  
Northampton NN2 6NE

Distributed by/Distribuido por/Distribué par/Vertrieb durch/Distribuito da:  
And in the USA/Y en EEUU/Et aux Etats-Unis/Und in USA/E in USA:

Denver, PA 17517 USA

Rev. 2013/12

**9793164**

Blue

Blu

Blue

Blau

Blu

Blauw

**9793165**

Grey

Gris

Gris

Gräu

Grigio

Grijs

**9793166**

White


Blanco

Blanc

Weiß

Bianco

Wit

## DIRECTIONS FOR USE for dentists/dental technicians

Use: For the manufacturing of the wax scaffolding of dental implants with plastic covered crowns, with the help of CAD/CAM system.

1. **General**

The forming waxes can be processed with different milling systems.

2. **Storage**

Protect from humidity and sunlight.

3. **Area of application:**

The forming wax is not designed for use in oral application. The material can be used for the preparation of dental crowns and bridges. Since it burns without residues, it is suitable for the preparation of molds in case of casting and ceramic molding processes.

4. **Processing**

The wax milled scaffoldings can be removed with the help of appropriate cutters or cutting discs.

5a. **Installation of the casting channel**

The setting and the selection of the diameter of the casting channel depends on the used devices and alloys or pressed ceramics.

5b. **Embedding**

Line the silicon or metal muffle with a liner. Mix and embed the embedding masse under vacuum observing the directions of use provided by the manufacturers of the embedding masses and alloys.

6. **Disposal**

The uncontaminated residues can be re-used for activities related to wax. Since the product is insoluble in water, it can be disposed together with normal household waste.

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